


# NiTAN Stone Basket

Our NiTAN high-quality stone retrieval baskets are designed for capturing and extracting urinary stones from the urinary tract and minimise the risk of injury to patients. Tipless stone baskets are produced using a special laser procedure and are perfect for use with a flexible ureteroscope to extract stones from the kidneys.



**Ergonomic UV-design handle for a secure hold**

**Super-elastic, excellent kink resistance and flexibility for greater scope**

**High radial force and perfect 1:1 transmission of rotational movements**

**Blue colour coding for excellent endoscopic visibility**

**Supplied in sterile blister pack, ready-for-use**



**NiTAN helical and tipless stone baskets are available in a range of designs.**

Order number	Size (FR/CH)	Wires	Working length cm	Basket Ø mm	Basket shape
48012456	2,5	4	90	12	helical
48012457	3	4	90	12	helical
48012458	3	4	120	12	helical
48012461	4,5	4	90	15	helical
48012462	4,5	5	90	15	helical
48012459	1,8	4	120	11	tipless
48012460	2,5	4	90	12	tipless